



BRASSELER
USA®

BY YOUR SIDE

 **EVERclean**
AIR LINE SETUP GUIDE



► Order online at BrasselerUSA.com or call 800.841.4522.
In Canada call 800.363.3838.



EVERCLEAN AIR LINE SET UP GUIDE

1



LOCATE AIR OUTLET & REMOVE NUT

Locate your **wall outlet for air**. The outlet could be under your sink or above the counter.

Remove the **nut** from the wall; this may require an adjustable wrench. The **nut** should contain a **gold metal ring**.

2



PREPARE THE FITTING

Using **Brass Poly Fitting** (Order No. 5028726U0), slide the **nut** onto the **cylinder end** of the brass fitting first, then slide on the **gold ring**.

3



INSERT THE FITTING

Insert the **Brass Poly Fitting** back into the **air outlet** cylindrical side down, with the black poly fitting facing up. Use an adjustable wrench to tighten the silver nut back onto the wall outlet.

PREPARE EVERCLEAN AIR LINE

Unscrew the **top brass nut** from the **Brass Poly Fitting** and slide it down the **EverClean Maintenance Unit's air line**. Flange side should **face down** towards the end of the air line.

INSERT AIR LINE

Insert the **air line** onto the **brass stem** of the **Brass Poly Fitting** as far as it will go, then tighten the **top brass nut** to the fitting. Hand tighten as much as possible, then hold the **bottom brass nut** in place with one adjustable wrench and tighten down the **top brass nut** with a second adjustable wrench. This will keep the fitting from spinning in place and prevents over tightening. Tighten down the **top brass nut** until there are only 1 or 2 visible threads between the fitting and the base. When you're done, the fitting should look like the circled image to the right.

TURN ON AIR LINE & TEST UNIT

Slowly turn the **air line valve** on (lefty loosey) and check for any air leaks. If found, check for proper connection and proper tightness. Next, turn on the **EverClean Maintenance Unit** and hold down the **Start Button** to enter the **Set Up** function. Cycle through settings and make sure **Compressed Air** is set to **on**. Cycle to **Air Pressure Settings** and check that the air pressure is set at or close to **4.5 bars**. Complete 2-3 **Test Cycles** to ensure the unit is working properly.

4



5



6



► Order online at BrasselerUSA.com or call **800.841.4522**.
In Canada call **800.363.3838**.